See discussions, stats, and author profiles for this publication at: https://www.researchgate.net/publication/231170482

Carbon disulfide charge exchange as a lowenergy ionization technique for hydrocarbon characterization

ARTICLE in ANALYTICAL CHEMISTRY · MARCH 1993

Impact Factor: 5.64 · DOI: 10.1021/ac00054a020

CITATIONS

28

READS

25

2 AUTHORS:



Chang S Hsu
Florida State University

151 PUBLICATIONS 1,672 CITATIONS

SEE PROFILE



Kuangnan Qian

ExxonMobil

69 PUBLICATIONS **2,506** CITATIONS

SEE PROFILE

CS₂ Charge Exchange as a Low-Energy Ionization Technique for Hydrocarbon Characterization

Chang S. Hsu* and Kuangnan Qian*

Exxon Research and Engineering Company, P.O. Box 998, Route 22 East, Annandale, New Jersey 08801

CS2*+ ions formed under low-energy electron ionization conditions readily undergo charge-exchange reactions with hydrocarbons having ionization potential less than 9.5 eV. CS₂ charge exchange (CS₂/CE) has many advantages over low-voltage electron-impact ionization (LV/EI) in the characterization of hydrocarbons. Small variations in electron beam energies do not affect the relative yields of the molecular ions. Higher sensitivities in CS2/CE over LV/EI facilitate highresolution accurate mass measurements. Uniform sensitivities over a wide range of compound types eliminate the need for calibration compounds for the quantitative determination of the components in a hydrocarbon mixture.

INTRODUCTION

Since its introduction by Field and Hastings in 1956, lowvoltage electron-impact ionization (LV/EI) has been commonly used in the petroleum industry for compound type analysis of aromatic hydrocarbons in fossil fuels. In this technique, the electron energy is generally set around 10 eV so that aromatic hydrocarbons, which have ionization potentials below 9.5 eV, can be selectively ionized with minimum fragmentation, leaving saturated hydrocarbons largely unionized. This technique, when combined with either high-2-5 or low-6-11 resolution mass spectrometry, provides quantitative measurements of molecular ions, which yield the compound type distribution as well as the carbon number distribution within each compound type (or homologous series) of a hydrocarbon mixture. We have recently applied LV/EI highresolution MS coupled with HPLC for the characterization of mixtures of heavy hydrocarbons. 12,13

There are, however, several disadvantages associated with the LV/EI technique. The sensitivity of the mass spectrometer is relatively low because of low electron density and small ionization cross sections under LV/EI conditions. The

† Current address: W. R. Grace & Co., Washington Research Center, 7379 Route 32, Columbia, MD 21044.

- Field, F. H.; Hastings, S. H. Anal. Chem. 1956, 28, 1248-1255.
 Johnson, B. H.; Aczel, T. Anal. Chem. 1967, 39, 682-685.
 Lumpkin, H. E.; Aczel, T. ACS Symp. Ser. 70 1978, 261-273.
- (4) Schmidt, C. E.; Sprecher, R. F.; Batts, B. D. Anal. Chem. 1987, 59,
- (5) Schmidt, C. E.; Sprecher, R. F. Novel Techniques in Fossil Fuel Mass Spectrometry; Ashe, T. R., Wood, K. V., Eds.; ASTM STP 1019,
- 1989; pp 116-132.(6) Lumpkin, H. E. Anal. Chem. 1958, 30, 321-325.
- (7) Kearns, G. L.; Maranowski, N. C.; Crable, G. F. Anal. Chem. 1959,
- (8) Crable, G. F.; Kearns, G. L.; Norris, M. S. Anal. Chem. 1960, 32,
- (9) Lumpkin, H. E.; Aczel, T. Anal. Chem. 1964, 36, 181-184.
- (10) Hsu, C. S.; Aczel, T.; Ashe, T. R.; Fedora, J. W. Proceedings of the 36th ASMS Conference on Mass Spectrometry and Allied Topics; San Franscisco, CA, June 5-10, 1988; pp 326-327.
 (11) Aczel, T.; Hsu, C. S. Int. J. Mass Spectrom. Ion Phys. 1989, 92,
- (12) Hsu, C. S.; McLean, M. A.; Qian, K.; Aczel, T; Blum, S. C.; Olmstead, W. N.; Kaplan, L. H.; Robbins, W. K.; Schulz, W. W. Energy Fuels 1991, 5, 395-398
 - (13) Qian, K.; Hsu, C. S. Anal. Chem. 1992, 64, 2327-2333.

effective electron energy is sensitive to ion source conditions, especially the space charge and cleanliness of the ion source. Small changes in electron beam energy can cause large variations in ion yield. In other words, in LV/EI conditions the ionization efficiency curve has a steep slope.

Other ionization techniques have also been used for hydrocarbon characterization. Field ionization and field desorption (FI/FD) have been gaining popularity in recent years 14-21 because of their unique features, such as uniform sensitivity across a broad range of compound types, applicability to both saturated and aromatic hydrocarbons, and the ability to measure less volatile hydrocarbons. FI has been used for volatile organics while FD has become the method of choice for the characterization of heavy hydrocarbons. However, with FI/FD the mass spectrometer can only be operated at low resolution and cannot at present be coupled with a liquid chromatograph for on-line liquid chromatography/mass spectrometric (LC/MS) analysis of heavy complex hydrocarbon mixtures. To characterize heavy hydrocarbons using FDMS, tedious preparative liquid chromatographic separations are usually required.

Chemical ionization (CI), such as proton transfer and charge exchange, has also been explored for hydrocarbon analysis. 22-27 Munson and co-workers²⁸ recently suggested the use of the benzene charge-exchange reaction for aromatic hydrocarbon analysis, which yields uniform molar sensitivity for various hydrocarbon compounds. This method was also applied to differentiate between paraffins and olefins: paraffins in general do not react with benzene ions by charge transfer whereas olefins do. However, the rather low volatility of benzene makes it difficult to introduce enough benzene vapor into a combined EI/CI source or a CI source for a moving belt interface, which has a finite opening to allow the belt moving into the ion source, to produce a stable stream of benzene ions, resulting in low sample ion yields. A specially designed

- (18) Boduszynski, M. M. Energy Fuels 1987, 1, 2-11. (19) Boduszynski, M. M. Energy Fuels 1988 2, 597-613.
- (20) Gellegos, E. J. Proceedings of the 37th ASMS Conference on Mass Spectrometry and Allied Topics; Miami Beach, FL, May 21-26, 1989; pp
- (21) Reynolds, S. D.; Aczel, T. Proceedings of the 33rd ASMS Conference on Mass Spectrometry and Allied Topics; San Diego, CA,
- May 26-31, 1985; pp 656-657. (22) Hatch, F.; Munson, B. Anal. Chem. 1977, 49, 731-733.
 - (23) Subba Rao, S. C.; Fenselau, C. Anal. Chem. 1978, 50, 511-515.
- (24) Hunt, D. F.; Harvey, T. M. Anal. Chem. 1975, 47, 2138-2141. (25) Fan, T. P.; Mango, F. D. Proceedings of the 39th ASMS Conference on Mass Spectrometry and Allied Topics; Nashville, TN, May 19-24, 1991; pp 1099-1100.
- (26) Wadsworth, P. A.; Villalanti, D. C. Hydrocarbon Process. 1992; May, 109-112.
- (27) Dzidic, I.; Peterson, H. A.; Wadsworth, P. A.; Hart, H. V. Anal. Chem. 1992, 64, 2227-2232 (28) Allgood, C.; Ma, Y. C.; Munson, B. Anal. Chem. 1991, 63, 721-725.

⁽¹⁴⁾ Scheppele, S. E.; Grizzle, P. L.; Greenwood, G. J.; Marriott, T. D.; Perreira, N. B. Anal. Chem. 1976, 48, 2105-2113.

⁽¹⁵⁾ Scheppele, S. E.; Hsu, C. S.; Marriott, T. D.; Benson, P. A.; Detwiler, K. N.; Perreira, N. B. Int. J. Mass Spectrom. Ion Phys. 1978, 28, 335-346. (16) Boduszynski, M. M.; Hurtubise, R. J.; Allen, T. W.; Silver, H. F. Anal. Chem. 1983, 55, 225-231.

⁽¹⁷⁾ Boduszynski, M. M.; Hurtubise, R. J.; Allen, T. W.; Silver, H. F. Anal. Chem. 1983, 55, 232-241.

or dedicated CI source with a heated reagent gas transfer line to effectively introduce benzene vapor into the ion source would be needed for the application of this technique.

Another reagent gas with an ionization potential similar to benzene's is CS₂, which has been used as a reagent gas for charge-exchange mass spectrometric studies.²⁹ CS₂ does not contain hydrogen so that proton transfer from CS2*+ is impossible. In addition, the CS2*+ ion is not a good electrophile to form an adduct ion with an organic molecule. Therefore, charge exchange is the predominant process for CS2*+ to interact with organic molecules with ionization potentials less than 9.5 eV, yielding simplified spectra of almost exclusively molecular ions. We have recently evaluated CS₂ charge exchange (CS₂/CE) as a low-energy ionization technique for hydrocarbon characterization. The fact that the boiling point of CS₂ is only 46 °C, 36 °C lower than that of benzene, facilitates charge-exchange experiments using a combined EI/CI source or a CI source for a moving belt interface without any modification.

EXPERIMENTAL SECTION

The comparison between CS_2 charge exchange and low-voltage EI ionization techniques was carried out using a VG-70 VSE double-focusing mass spectrometer which is coupled with a Hewlett-Packard 5890 gas chromatograph (GC) and a Varian 5560 ternary liquid chromatograph (LC) via a moving belt interface. GC/MS was used for the studies of relative sensitivities of polynuclear aromatic hydrocarbons and alkylbenzenes. In GC/MS, a 25-m \times 0.25-mm-internal diameter fused silica capillary column coated with 0.25- μ m HP-1 methyl silicone stationary phase was used. The column temperature was programmed to rise from 100 to 300 °C at 10 °C/min. Helium was used as a carrier gas.

In the LC/MS experiments, a petroleum fraction boiling between 650 and 1050 °F was selected to compare the sensitivities of CS $_2$ /CE and LV/EI ionization techniques. The LC/MS conditions have been described previously. 12 In brief, the petroleum fraction was separated by HPLC into elution regions of saturates, monoaromatics, diaromatics, triaromatics, tetraaromatics, and polars. The LC eluants were continuously transported into the ion source of the VG-70 VSE double-focusing mass spectrometer using a moving belt interface. The resolution of the mass spectrometer was set at 10 000. The nose heater current of the moving belt interface was set at 2.4 A with an estimated temperature of 300 °C.

To set up LV/EI conditions, xylene was used as a reference compound to set the electron beam energy and tune the mass spectrometer. The electron beam voltage was adjusted until the ratio of the intensities of the m/z (mass-to-charge) 106 molecular ion to the m/z 91 fragment ion of xylene was greater than 50:1. This ratio matches the largest achievable ratio in CS₂/CE. For CS₂/CE experiments, CS₂ was introduced into the ion source of the VG-70 VSE mass spectrometer via a standard chemical ionization (CI) reagent gas inlet. Since there was no pressure measuring device in the ion source, the source pressure is adjusted until the source chamber gauge read 2×10^{-5} Torr. The actual pressure inside the ion source is estimated to be in the range 0.1–0.5 Torr.

 CS_2 was purchased from EM Science with a purity of 99.9%. Xylene from J. T. Baker Chemical Co. Polynuclear aromatic hydrocarbons and alkylbenzenes were purchased from Aldrich.

RESULTS AND DISCUSSION

In low-voltage electron-impact ionization, molecular ions are formed with a distribution of internal energies. In contrast, in CS_2 charge exchange the internal energy of the molecular ion is well defined and equal to the difference between the recombination energy (RE) of $CS_2^{\bullet+}$ cation and the ionization potential (IP) of the sample molecule.³⁰

$$CS_2^{\bullet +} + M \rightarrow CS_2 + M^{\bullet +}$$

$$\Delta E = RE(CS_2^{\bullet +}) - IP(M)$$

However, almost all of the CS2 charge-exchange experiments reported in the literature used CS2*+ ions formed by highenergy ionization processes, such as by reaction with N2*+ or by direct bombardment with 70-eV electrons.²⁹ Under these high-energy (and high pressure for chemical ionization) conditions, fragment and cluster ions, such as S_2 *+ and (CS₂)S+ with different recombination energies, are also produced. CS₂*+ itself may also be electronically and/or vibrationally excited. Thus, the charge-exchange process is not limited to $CS_2^{\bullet+}$ ions in the ground state. Hence, to exclude these other ionic species and higher energy states of CS2*+, it is desirable to ionize CS₂ using a low-energy electron beam. Table I lists the relative abundances of these species for electron beam energies of 70 and 10 eV. The multiple ionic species formed at 70 eV can result in complicated charge-exchange spectra. At 10 eV, essentially all of the ionic species formed are CS₂⁺

To establish optimal LV/EI conditions, xylene is commonly used as a reference compound. To illustrate the dependence of the degree of fragmentation on electron energy, the relative abundances of the molecular ion (m/z 106) and two fragment ions (m/z 91 and 15) of xylene are plotted versus nominal electron beam energy in Figure 1. To minimize fragmentation the nominal electron energy is usually set below 10 eV, which is accompanied by a significant loss in ion yield, as shown in Figure 2. Chemical ionization (CI) generally has higher sensitivity than electron-impact ionization (EI) considering the rate constant in CI and ionization cross sections in EI.²⁹ For this reason CS₂/CE is expected to provide higher sensitivity than LV/EI; this expectation is borne out in Figure 2 over an electron energy range of 9-16 eV. Specifically, in the range of 9-10 eV which is commonly used in LV/EI, CS₂/ CE is about 50 times more sensitive.

Because of this improvement in MS sensitivity (i.e., signal-to-noise ratio) we realize an increase in mass measurement precision. An example is shown in Table II which lists measured mass defects of four ions, $C_{33}H_{50}$ (MW 446.391), $C_{33}H_{52}$ (MW 448.407), $C_{31}H_{46}S$ (MW 450.332), and $C_{31}H_{48}S$ (MW 452.348), by CS_2/CE and LV/EI. Averaged over 10 consecutive scans, both techniques yield similar results in mass measurement accuracy (1–2 millimass units (mmu) from the actual masses). However, the increased sensitivity of the mass measurement in CS_2/CE allows this technique to be more precise than LV/EI with the same number of scans: 1 versus 3 mmu. With this high precision fewer scans are needed in CS_2/CE for accurate mass measurement, which is an important feature for components in low abundance (weak peaks) and for samples with limited quantities.

The dependence of the peak intensity ratio, m/z 91 to m/z 106, on nominal electron beam energy in CS_2/CE and LV/EI is illustrated in Figure 3. In LV/EI this ratio increases rapidly from 1% to 65% as the electron beam energy increases from 9 to 16 eV. In contrast, it increases very slowly from 1% and maximizes at around 7% in CS_2/CE .

To achieve a desirable condition for either LV/EI or CS₂/CE, i.e., maximum sensitivity with minimum fragmentation, the electron beam energy is reduced to around 10 eV so that the ratio of the m/z 106 ion to the m/z 91 ion is greater than 50:1. This ratio is much less sensitive to the electron beam energy in CS₂/CE than in LV/EI, thereby greatly simplifying the tuning procedure and improving reproducibility.

Table I. Mass Spectral Data of CS2 Charge-Exchange Gas Ionized by an Electron Beam at 70 and 10 eV

	ionic species	fractional abundance		
m/z^a		70 eV	10 eV	
32	S+	0.83		
44	CS ⁺	0.54		
64	$\mathbf{S_2}^+$	11.64		
76	$\overset{\circ}{\mathbf{CS}_2}^+$ $\overset{\circ}{\mathbf{S}_3}^+$	79.34	99.01	
96	\mathbf{S}_3 +	0.89		
108	$(CS_2)S^+$	2.03		
152	$(CS_2)_2^+$	1.06	0.99	
160	S ₅ +	2.82		
192	$(CS_2)_2^+ S_5^+ S_6^+$	0.85		
sum	· ·	100	100	

a m/z of the lightest isotopic peak.




Figure 1. Relative abundances of xylene molecular and fragment ions as a function of electron beam energy in LV/EI.




Figure 2. Ion yield of the xylene molecular ion (m/z) 106) as a function of electron beam energy in LV/EI and CS₂/CE.

One of the difficulties of using LV/EI for the quantification of components in hydrocarbon mixtures is the scarcity of alkylated aromatic and hydroaromatic model compounds for sensitivity calibration. This calibration is important because the ionization cross section in EI strongly depends on the polarity and structure of the molecules. This difficulty can be alleviated to some extent by using CS₂/CE ionization. The relative molar sensitivities of alkylbenzenes (a compound series with various lengths of alkyl substituents) and polynuclear aromatic hydrocarbons (PAH's, core structures with different degrees of aromaticity) in CS2/CE and LV/EI, at a source temperature of 280 °C, are shown in Figures 4 and 5. Clearly, the relative sensitivity in CS2/CE is largely independent of molecular structure as in benzene charge-exchange ionization.²⁸ In LV/EI, on the other hand, relative molar sensitivities are strongly dependent on molecular weight and molecular structure. 4,6,7

It is often reported in the literature that the temperature of the ion source be kept below 200 °C even for quantitative analysis of rather high-boiling hydrocarbon mixtures. 1-23,24,27,28 This could have been due to the consideration of fragmen-

Table II. Mass Defects of Four Ions at Nominal Masses of 446, 448, 450, and 452 Measured by CS2/CE and LV/EI

		nominal mass				
	scan no.	446	448	450	452	
CS ₂ /CE	1	0.392	0.406	0.339	0.348	
	$\frac{2}{3}$	0.390	0.406	0.334	0.348	
		0.392	0.405	0.334	0.348	
	4	0.391	0.405	0.333	0.348	
	5	0.392	0.408	0.334	0.347	
	6	0.393	0.405	0.334	0.350	
	7	0.392	0.405	0.333	0.348	
	8	0.395	0.408	0.334	0.350	
	9	0.392	0.408	0.332	0.349	
	10	0.394	0.407	0.335	0.350	
	mean	0.392	0.406	0.334	0.349	
	std dev	0.001	0.001	0.002	0.001	
LV/EI	1	0.390	0.402	0.335	0.345	
	2	0.393	0.407	0.333	0.349	
	3	0.392	0.408	0.332	0.350	
	4	0.392	0.409	0.335	0.347	
	5	0.388	0.407	0.330	0.343	
	6	0.390	0.404	0.326	0.346	
	7	0.394	0.403	0.335	0.349	
	8	0.397	0.403	0.334	0.349	
	9	0.398	0.407	0.337	0.346	
	10	0.396	0.406	0.335	0.344	
	mean	0.393	0.408	0.333	0.347	
	std dev	0.003	0.004	0.003	0.002	
actual mass defect		0.391	0.407	0.332	0.348	




Figure 3. The intensity ratio of the m/z 91 to m/z 106 ion of xylene as a function of electron beam energy in LV/EI and CS $_2$ /CE.




Figure 4. Relative molar sensitivity of alkylbenzenes as a function of molecular weight in LV/EI and CS₂/CE.

tation of molecular ions caused by excess thermal energy at higher source temperatures. Figure 6 shows the dependence of the relative sensitivity of PAH's in CS₂/CE on the source temperature. At a source temperature of 280 °C the molar sensitivities of all PAH's are almost equal. However, the sensitivities of heavy PAH's decrease rapidly as the source block temperature decreases from 280 to 150 °C. Similar phenomena were also observed in other ionization techniques; thus, one may conclude that this loss in sensitivity for heavier




Figure 5. Relative molar sensitivity of polynuclear aromatic hydrocarbons as a function of molecular weight in LV/EI and CS₂/CE.




Figure 6. Relative molar sensitivity of polynuclear aromatic hydrocarbons in CS₂/CE as a function of ion source temperature.

PAH's is mainly due to their condensation in a cold source. For a good result in quantification, therefore, a preferred ion source temperature would be above 280 °C.

LC can be easily combined with LV/EI and CS₂/CE MS via a moving belt interface. The characterization of heavy hydrocarbons by LC/MS under LV/EI conditions has been discussed previously.^{12,13} The relative amounts of sample in the various LC elution regions as measured by HPLC using an evaporative mass detector (EMD), LC/MS with LV/EI and LC/MS with CS₂/CE, are shown in Figure 7. The results by EMD is expressed as weight % after sensitivity correction, whereas those by LC/MS are percentages of total ionization without sensitivity corrections. CS2/CE shows much better agreement with EMD than does LV/EI because of its uniform sensitivities for all aromatic compound types. Using EMD results as a reference, the tendency of overestimating 2+ring aromatics and underestimating 1-ring aromatics by LV/ EI is evident. This result is consistent with the fact that LV/EI sensitivity increases with the number of aromatic rings as shown in Figure 5.

Despite the rather large differences in overall sensitivities, the distributions of compound types within each LC elution region as determined by LV/EI and CS₂/CE show only small differences. The compound type distributions of monoaromatic hydrocarbons and diaromatic thiophenes of a petroleum distillate as determined by LV/EI and CS₂/CE are shown in Figures 8 and 9. In Figure 8 the C_nH_{2n-6} series are alkylbenzenes, the C_nH_{2n-8} series are mononaphthenobenzenes, the C_nH_{2n-10} are dinaphthenobenzenes, etc. It can be observed in this figure that LV/EI relative sensitivities of hydroaromatic hydrocarbons (C_nH_{2n+z} with z < -6) decrease with the number of naphthene rings attached to the aromatic core. For benzothiophenes ($C_nH_{2n-10}S$) and naphthenobenzothiophenes ($C_nH_{2n+z}S$ with z < -10) shown in Figure 9, the relative sensitivities are almost identical regardless of the number of naphthene rings.




Figure 7. Distribution of aromatic hydrocarbons determined by EMD, CS₂/CE, and LV/EI.




Figure 8. Distribution of hydrocarbon types in monoaromatic LC region determined by LV/EI and CS₂/CE.





Figure 9. Distribution of sulfur compound types in diaromatic LC region determined by LV/EI and CS $_2$ /CE.

The carbon number distributions of each compound series determined by LV/EI and CS_2/CE in general show very good agreement. An example is shown in Figure 10 for dinapthenobenzenes (C_nH_{2n-10} series in monoaromatic LC elution region). Both techniques yield a carbon number distribution centered around C_{34} with a full width at half height of about 12 carbons.

CONCLUSIONS

Charge exchange using CS_2 as a reagent gas (CS_2/CE) has been evaluated and compared with a traditional low-voltage electron-impact ionization (LV/EI) technique for the characterization of hydrocarbons. It would be desirable to obtain primary $CS_2^{\bullet+}$ ion using a low-energy electron beam. The optimal ionization conditions for CS_2/CE are much easier to obtain than LV/EI because ionization of the molecules depend on the well-defined recombination energy of the $CS_2^{\bullet+}$ ions rather than the energy settings of the electron beam. Thus, small variations in electron beam energy do not affect the relative sensitivity significantly. In addition, CS_2/CE is found






Figure 10. Carbon number distribution of dinaphthenobenzenes determined by LV/EI and CS2/CE.

to be much more sensitive than LV/EI and therefore it provides a better mass measurement precision than LV/EI. This feature is particular useful when high resolution mass spectrometry is used for accurate mass measurements. Furthermore, CS₂/CE also provides more uniform molar sensitivity for the quantitative measurements of hydrocarbon molecules with different aromatic structures than LV/EI. Thus, it eliminates the necessities of obtaining a wide range of alkylated aromatic and hydroaromatic hydrocarbons for the calibration of their relative molar sensitivities for quan-

titation. Therefore, for the characterization of complex hydrocarbon mixtures CS₂/CE is recommended over LV/EI.

ACKNOWLEDGMENT

The valuable comments from Peter B. Grosshans during manuscript preparation are highly appreciated.

RECEIVED for review July 16, 1992. Accepted November 25, 1992.